



D-10[™]

for HbA_{1c}, HbA₂ and HbF





D-10™

for HbA_{1c}, HbA₂ and HbF


Flexible
to solve more testing needs

Comprehensive
β-thalassemia and
diabetes testing


Easy
for simple operation

Dependability is about more than keeping your laboratory running smoothly; it's about the quality diabetes care you support. That's why we developed the D-10™ with reliability and efficiency in mind.

A simple, fully-automated solution, the D-10™ combines diabetes and β-thalassemia testing, enabling rapid HbA_{1c} or HbA₂/F/A_{1c} testing using primary tube sampling—so you can accomplish more in fewer steps. With the D-10™, it's easier to deliver a full picture of diabetes treatment progress—and that can be the difference for the people who count on you most.



Diabetic patient: HbA_{1c} program



HbS: 31%, a presumptive HbAS patient
HbA₂: Quantified at 3.1%
HbA_{1c}: Reportable in the presence of HbS

Experience a new standard in flexibility

- Same reagent kit for:
- HbA_{1c} Short Program
 - Dual Program for HbA_{1c} as well as screening for abnormal hemoglobins
 - Quantitation of HbA₂ and HbF levels aid in the diagnosis for β-thalassemia

quality results



D-10™

for HbA_{1c}, HbA₂ and HbF

Flexible
to solve more testing needs

Comprehensive
β-thalassemia and
diabetes testing

Easy
for simple operation

Dependability is about more than keeping your laboratory running smoothly; it's about the quality diabetes care you support. That's why we developed the D-10™ with reliability and efficiency in mind.

A simple, fully-automated solution, the D-10™ combines diabetes and β-thalassemia testing, enabling rapid HbA_{1c} or HbA₂/F/A_{1c} testing using primary tube sampling—so you can accomplish more in fewer steps. With the D-10™, it's easier to deliver a full picture of diabetes treatment progress—and that can be the difference for the people who count on you most.



Where effective diabetes and thalassemia treatment begins

Because completing multiple tests in fewer steps supports more than just your workflow—it's a big part of delivering timely, effective care.

Flexible


- Easy switching between HbA_{1c} and β-thalassemia testing without changing reagents or cartridge
- Rapid 3 minute HbA_{1c} program
- Complete β-thalassemia program including HbA_{1c} results
- Expandable for higher volumes with the D-10™ Rack Loader

Comprehensive Testing

- β-thalassemia and diabetes testing
- Simultaneous quantitative HbA₂, HbF and HbA_{1c} determination
- Accurate HbA₂ and HbF results in just 6.5 minutes
- Detection of common hemoglobin variants
- HbA_{1c} reportability in the presence of Hb S, C, D and E

Easy to Use

- Touchscreen operation
- One-step cartridge switching
- Fully automated with barcode reading and primary tube sampling
- Simplified user maintenance



Powerful results and the flexibility you need to be the difference

Like you, we take quality patient care seriously. We developed the D-10™ to deliver complete and accurate hemoglobin testing in minutes. With the flexibility to run a HbA_{1c} or a HbA₂/F/A_{1c} test on one platform with only one set of reagents and only one cartridge, the D-10™ allows laboratories to identify patients with β-thalassemia or abnormal hemoglobins. A comprehensive, easy solution, the D-10™ can help you be the difference in patient care.

*Every patient deserves quality results. You deserve **flexibility**.*

D-10™ for HbA_{1c}, HbA₂ and HbF

Unparalleled Support

Bio-Rad's world class service and support team is backed by our decades-long leadership in the field of HbA_{1c} testing. The D-10™ represents our continued commitment to helping you be the difference for the patients who count on you.

Powerful testing, simplified

D-10™ for HbA_{1c} or HbA₂/F/A_{1c}

Accomplish more in fewer steps with the D-10™'s simple operation, easy switching between HbA_{1c} and β-thalassemia assays and expandability for higher volumes.



Ordering Information

Catalog No.	Description
220-0220	D-10™ Hemoglobin Testing System (110/220V) with Accessories.....1 unit
220-0600	D-10™ Rack Loader.....1 unit
220-0101	D-10™ Hemoglobin A _{1c} Program Reorder Pack.....400 tests
220-0201	D-10™ Dual Program Reorder Pack400 HbA _{1c} tests or 200 HbA ₂ /F/A _{1c} tests
220-0375	D-10™ Printer Paper10 rolls
740	Lyphochek® Diabetes Control Bilevel (3 each of 2 levels).....6 x 0.5 mL
171	Liquichek™ Diabetes Control Level 16 x 1 mL
172	Liquichek™ Diabetes Control Level 26 x 1 mL
173	Liquichek™ Diabetes Control Level 36 x 1 mL
553	Lyphochek® Hemoglobin A ₂ Control Bilevel (2 each of 2 levels)4 x 1 mL
120	Lyphochek® Hemoglobin A _{1c} Linearity Set (1 each of 4 levels).....4 x 0.5 mL
BC80	EQAS Hemoglobin Program (3 analytes).....12 x 0.5 mL

Bio-Rad **A1c** • Be the difference



**Bio-Rad
Laboratories**

For further information, please contact the Bio-Rad office nearest you or visit our website at www.bio-rad.com/diagnostics

Clinical
Diagnostics Group

Website www.bio-rad.com/diagnostics U.S. 1-800-2BIO-RAD Australia 61-2-9914-2800 Austria 43-1-877-8901 Belgium 32-9-385-5511 Brazil 5521-3237-9400
Canada 1-514-334-4372 China 86-21-64260808 Czech Republic 420-241-430-532 Denmark +45-4452-1000 Finland 358-9-804-22-00 France 33-1-47-95-60-00
Germany +49-(0)89-318-840 Greece 30-210-7774396 Hong Kong 852-2789-3300 Hungary 36-1-455-8800 India 91-124-4029300 Israel 972-3-9636050 Italy
+39-02-216091 Japan 81-3-6361-7070 Korea 82-2-3473-4460 Mexico +52(55)5488-7670 The Netherlands +31-318-540666 New Zealand 64-9-415-2280 Norway
47-23-38-41-30 Poland 48-22-3319999 Portugal 351-21-472-7700 Russia 7-495-721-14-04 Singapore 65-6415-3188 South Africa 27-11-442-85-08 Spain
34-91-590-5200 Sweden 46-8-555-127-00 Switzerland 41-61-717-95-55 Thailand 662-651-8311 United Kingdom +44-(0)20-8328-2000