

Easy Wash

Microplate Washer

Control

Fully Automatic

Accuracy

Reliable Result

Technology


Multi-assay feature

Performance

Ability to communicate





Features



- ✓ Fully automatic.
- ✓ Easy Windows operation system, large LCD display.
- ✓ Washes flat, round and V-bottom strips and plates.
- ✓ 12-way or 8-way manifold, automatic and manual positioning.
- ✓ Large memory to store up to 50 wash protocols.
- ✓ Automatic monitoring for vacuum and pressure, automatic rinse cycle.
- ✓ Automatic check for plate and strip.
- ✓ Washing pressure, volume and time adjustable.
- ✓ Liquid level sensor and alert feature.
- ✓ No special airtight bottle required for wash bottle.
- ✓ Easy maintenance design.

Easy Wash


Microplate Washer




Main Specifications

Plate Mode	96-or 48- well plate or strip
Dispense Precision	CV<2%
Residual Volume	< 1ul
Washing Volume	50-3000ul for 8-way manifold 50-2000ul for 12-way manifold
Washing Cycle	Maximum 99 cycles
Strip Washing	1 to 12 channels adjustable
Soaking Time	0-24 hr
Display	5" LCD
Bottle Capacity	2 litres each
Memory	50 washing protocols
Interface	RS-232 serial port
Net Weight	8KG
Dimensions LXWXH (mm)	390 x 340 x 180
	AC 220V \pm 10%, 50-60Hz, 110V on request


Program interface



Select strip to wash



Manifold positioning



System setting

