

GM Test



FungiXpert[®] **Aspergillus Galactomannan ELISA Detection Kit**

Product code: FGM096-001

Obtain result in 2h

Product Advantages

More Advance

International leading edge detection method, high sensitivity and specificity

More Accurate

Optimize the operation step, reduce the risk of contamination during the experiment

Faster

One-step detection, reducing the number of incubation and washing

Obtain results within 2 hours

More Economical

Split microplate, cost saving

Clinical Performance

- **Low detection limit**

Limit of Detection (LoD) is 0.5 ng/ml

- **Sensitivity and specificity**

	Serum	BAL
Sensitivity	89.8%	92.4%
Specificity	90.7%	94.1%

Product Characteristics

Basic Parameters	Aspergillus Galactomannan ELISA Detection Kit
Method	ELISA
Sample Type	Serum, BAL fluid
Sample Volume	300 µL
Specification	96 tests/kit
Detection Time	2 h
Detection Objects	<i>Aspergillus spp.</i>
Stability	The kit is stable for 1 years at 2°C-8°C
Limit of Detection	0.5 ng/ml

- **Recommendations**

Recommended by IDSA guideline for Aspergillosis 2016 (Patterson et al., 2016) and ESCMID-ECMM-ERS guideline for Aspergillosis 2018 (Ullmann et al., 2017)

For more information contact :



16, G. F. Arjun Arcade Phase - 1, Shopping Complex, Sector 30 & 33
Faridabad - 121003, Haryana (Delhi/NCR)
Mob. : +91 9654454435, 9654454425
Email. : nitin.gupta@bio-state.com Website. : www.bio-state.com